


Functional microbiome research – bioinformatics section

Part 2 – Metagenomics


You can grab today's scripts here:

```
wget https://raw.githubusercontent.com/danchurch/FunctionalMicrobiomePractical2022/main/funmic2023/funMetagenomicScript.txt
```




This morning we will attempt to:

- look at the quality of our illumina and nanopore reads
- profile the members of the microbial community
- assemble our illumina-based metagenome
- check the quality of our nanopore and illumina assemblies




This morning we will attempt to:


- look at the quality of our illumina and nanopore reads
- profile the members of the microbial community
- assemble our illumina-based metagenome
- check the quality of our nanopore and illumina assemblies




Check our reads with FASTQC.





Illumina data often has two, paired-end read files. Why?



Let's look at what FASTQC can tell us.




Illumina



Nanopore

Metagenome assembly: we use two assemblers:

- MetaFLYE can handle the longer reads, shallower depth and higher error rates of Nanopore data.
- MEGAHIT is optimized for short reads with high accuracy (Illumina data).



MetaPhlan and biobakery – marker-based discovery of who is in our metagenome? Not just 16s!! MetaPhlan is still only available for illumina data.

